

CHU | UVC
BRUGMANN


IMAGING PORTAL

Short manual



[Information](#)

CHU|UVC BRUGMANN



Accessing all records of a visit.

Enter the Access Code and Control Code.
This information can be found on the document you
received during your visit at the medical imaging
department.

Access Code

Control Code

ENTER

HOW TO CONNECT

Connect to the Imaging Portal Site : <https://imaging.chu-brugmann.be> and enter the access code and control code mentioned on the document you received during your visit to Imaging Department, or on the exam report.

CHU|UVC BRUGMANN

Information

CHU|UVC BRUGMANN

Accessing all records of a visit.

Enter the Access Code and Control Code.
This information can be found on the document you received during your visit at the medical imaging department.

Access Code

Control Code


ENTER

AGFA HealthCare

Browser Compatibility List

Copyright (c) 2015 Agfa HealthCare N.V., Sepestraat 27, B-2040, Morselt, Belgium.

MAIN MENU



STUDY OPTIONS

Search | PODNIK, PETE | 27 Male 01-Jan-1933 | Help | Logout

CT LUMBAL WZ
09-Oct-2001 11:56 (63065)

2: 3D

REPORTS

- Study Comments ← Display reports

LAYOUT

- 1x1 Layout
- 2x1 Layout
- 2x2 Layout

SHARE STUDY

- Share Study ← Obtain a QR code to share with a mobile device

IMPORTANT IMAGES

- For Patient ← Show key images only

Home

105 of 249

IMAGE SERIES NAVIGATION

Search PODNIK, PETE | 27 Male 01-Jan-1933 Help | Logout

CT LUMBAL WZ
09-Oct-2001 11:56 (63065)

2: 3D

Number of Series

Home To display a series, tap or drag a series from the series tray

Currently displayed series

1: Topogram
2: 3D
3: SAG
4: SAG (Series 3)

105 of 249

-646.00

The screenshot shows a medical imaging software interface. At the top, there is a search bar with the text 'PODNIK, PETE | 27 Male 01-Jan-1933' and a 'Help | Logout' link. Below this, the study name 'CT LUMBAL WZ' and date/time '09-Oct-2001 11:56 (63065)' are displayed. The main area shows a 3D axial CT scan of the lumbar spine. A 'Home' button is overlaid on the scan. To the right, a 'series tray' lists four series: '1: Topogram', '2: 3D', '3: SAG', and '4: SAG (Series 3)'. An arrow points from the '4' icon in the top right corner to the 'Number of Series' text. Another arrow points from the '4: SAG (Series 3)' item in the series tray to the 'Currently displayed series' text. The bottom of the interface shows a navigation bar with '105 of 249' and a value '-646.00'.

IMAGE NAVIGATION

Search | PODNIK, PETE | 27 Male 01-Jan-1933 | Help | Logout

CT LUMBAL WZ
09-Oct-2001 11:56 (63065)

2: 3D

Home

R L

Navigate back one image

Drag to navigate series

Image Number

Navigate to next image

Play

105 of 249

-646.00

The screenshot displays a medical software interface for image navigation. At the top, a search bar shows the patient name 'PODNIK, PETE' and date of birth '27 Male 01-Jan-1933'. Below this, the study name 'CT LUMBAL WZ' and acquisition details '09-Oct-2001 11:56 (63065)' are visible. The main window shows a 2D axial CT scan of the lumbar spine. A 'Home' button is overlaid on the image. On the left, a vertical toolbar contains icons for pan, zoom, and other navigation functions. At the bottom, a navigation bar includes a 'Play' button, a left arrow for 'Navigate back one image', a slider for 'Drag to navigate series' (currently at '105 of 249'), and a right arrow for 'Navigate to next image'. A numerical value '-646.00' is shown in the bottom right corner of the navigation bar.

CINEVIEW MODE NAVIGATION

The screenshot displays a medical software interface for viewing a CT scan in cineview mode. The main window shows a cross-sectional CT scan of the lumbar spine. The interface includes several control elements:

- Top Bar:** Search field, patient information (PODNIK, PETE | 27 Male 01-Jan-1933), Help, and Logout buttons.
- CT Scan Header:** CT LUMBAL WZ, 09-Oct-2001 11:56 (63065).
- Left Panel:** 2: 3D view, navigation icons (pan, zoom, reset), and a 'Rock Loop' control.
- Center:** The CT scan image with a 'Home' button overlaid.
- Right Panel:** Playback speed adjustment slider and a 'Slow frame rate warning' icon.
- Bottom Panel:** Playback controls including a 'Pause' button, 'Left trim guide', 'Drag to navigate series' (a blue bar with a slider), 'Right trim guide', and a '105 of 249' frame counter.

Annotations with arrows point to the following features:

- Rock
- Loop
- Pause
- Left trim guide
- Drag to navigate series
- Right trim guide
- Slow frame rate warning
- Playback speed adjustment

IMAGE TOOLSET

The screenshot displays a medical imaging software interface. At the top, a search bar contains the text "PODNIK, PETE | 27 Male 01-Jan-1933". Below this, the patient information "CT LUMBAL WZ" and "09-Oct-2001 11:56 (63065)" is visible. The main area shows a 3D axial CT scan of the lumbar spine. On the left side, there is a vertical toolbar with several icons. White arrows point from text labels to these icons. The labels are: "Flip or rotate images" (pointing to a four-way arrow icon), "Adjust Brightness/Contrast (Window Level)" (pointing to a sun icon), "Zoom or pan" (pointing to a magnifying glass icon), "Caliper measurements" (pointing to a caliper icon), "Link series" (pointing to a circular arrow icon), and "Reset all transformations" (pointing to a square with a circular arrow icon). At the bottom of the interface, there is a "Home" button and a page indicator "105 of 249".

Search PODNIK, PETE | 27 Male 01-Jan-1933 Help | Logout

CT LUMBAL WZ
09-Oct-2001 11:56 (63065)

2: 3D

Flip or rotate images

Adjust Brightness/Contrast (Window Level)

Zoom or pan

Caliper measurements

Link series

Reset all transformations


Home

105 of 249

PATIENT ORIENTATION INDICATORS



A = anterior
P = posterior
H = head
F = foot
L = left
R = right



A complete help section and some video tutorials are available by clicking the « Help » button in the upper right corner of the image viewer.

Not all the features mentioned on the help pages are available in the Imaging Portal of the C.H.U. Brugmann.